單元0:

虚擬儀控簡介

主題:

- a. LabVIEW的歷史
- b. 虛擬儀控
- c. LabVIEW優點列舉

LabVIEW的歷史



LabVIEW 8.0 Project Window

May 2003

LabVIEW 7 Express VIs, I/O Assistants, FPGA/PDA targets

August 2000

LabVIEW 6/Internet-ready measurement intelligence

March 1998

LabVIEW 5.0 ActiveX, Multithreading

February 1996

LabVIEW 4.0 Added professional tools, improved debugging

August 1993

LabVIEW 3.0 Multiplatform version of LabVIEW

September 1992

LabVIEW for Windows

January 1990

LabVIEW 2.0 for Macintosh


October 1986

LabVIEW 1.0 for Macintosh

April 1983

LabVIEW project begins

工業界標準的開發工具



Source: Survey of 400 US readers from T&M World, EDN, Design News, and R&D magazines, Q1 2004

支援各式各樣的軟硬體

其他軟體

- Wolfram Research Mathematica®
- ■Microsoft Excel®
- MathSoft MathCAD®
- Electronic Workbench MultiSim®
- Texas Instruments Code Composer Studio®
- The MathWorks MATLAB® and Simulink®
- ■Ansoft RF circuit design 軟體
- Microsoft Access®
- Microsoft SQL Server®
- Oracle[®]







通訊協定

- Ethernet
- ■CAN
- TCP/IP
- DeviceNet
- ■USB

- IFFF 1394
- **RS-232**
- **GPIB**
- RS-485
- 藍芽
- 其他

支援的作業系統

- Microsoft Windows
- Apple Macintosh
- Linux







利用普及化的商用電子技術

軟體技術 LabVIEW ※ 然 Windows 作業 MAC Java 程式語言 Basic 2000s 1980s 1990s




Virtual Instrumentation Defined


Virtual instrument system (n.):

Software-defined system, where software based on user requirements defines the functionality of generic measurement hardware


Virtual Instrumentation Explained



虚擬儀控組成



虚擬儀控架構



量測與自動化軟體架構

測試及資料管理軟體

NI TestStand , DIAdem

互動式工具

SignalExpress VI Logger Vision Builder Al 應用開發環境

LabVIEW

圖形化開發

量測自動化服務軟體





Signal Conditioning















所有的解決方案使用同一個平台




模組化開發環境

- ■模組化的設計
- ■利用副程式建立階層式的系統
- ■可重覆使用的程式碼



簡化電腦語言架構



利用圖控式語言快速完成開發


LabVIEW * Copyright (c) 1993, 1994 Ron Winacott * This program may be used, copied, modified, and redistributed freely * for ANY purpose, so long as this notice remains intact. */ ıı 🕐 ms ❿ Ōw 1 kHz dt [1000] t0 0 d = optarg; fprintf(stderr,"tgrep: Debug on at level(s) ") while (*d) { 3₂₁ 100 Waveform Graph DAQ Assistant L7660 data M.D.J.S flace & FE REGSEP;
{
 yet oction !;
 juint(steel "times East supply a seed in string(s) '
 warge()' 15 (Stop pool nov */ fl-q (flags & FE PEGSEP) {
 jorgs octude ();
 fpronts(syders "tgres, Kust supply a second cold ();
 sud file list is directory's),
 ucsgo(); 45 (6) 25 16 ns 🕐 🕾 **①** $\begin{array}{ll} & \text{fprintf}(\text{anderr} \ ^2\text{tgho}_{\Gamma} - \text{Kust supply a search atring}(\pi)) \\ & \text{sead file list ar directory}(\pi)). \end{array}$ Ōw 1 kHz poka dt [1000 norus 11 / 4.5 t0 0 iR af Waveform Graph fl-a 3₂, 100 **DAQ Assistant** ncpu: G. L7660 istatus = DAQ data if (tlade's SE REMEXE) {
 tor (s=0: s/regemp_tm; s++)
 pm_pot[i] = mokepot(string[s] NULE)
 use pmstch = 1; M.D. ulCount, piB iRetVal = NIDAOError ner 1: else (
____bet = bs_makepat(string[c]), -* only one ellowed *iIgnoreWarning); firms - FE COUNT iStatus = DAQ_Clear Stop iRetVal = NIDAQPlotWaveform printf(" The data is avai" search the limit - new .pen files - BOLD FDS - delaw file.
cascade_thr_thm: * search_thr_thm: * ,
* the number of files that can be open */
correct open files - search thr limit pshreed_astr_anit(&dctoshed_stsr).

LabVIEW優點列舉(1/6):相容DLL檔

```
(1) e:\...\séminaire RT\démonstrations\démo CVI\madII.c
File Edit View Build Run Instrument Library Tools Window Options Help
                                                                                 🔀 madll. vi Diagram
                                                                                                                                                   _ 🗆 ×
Operate Tools Browse Window Help
  #include <dataacg.h>
  #include <cvirte.h>
  static double mesure;
                                                                                                                           13pt Application Font
  int __stdcall DllMain (HINSTANCE hinstDLL, DWORD fdwReason, LPVO
      switch (fdwReason)
         case DLL_PROCESS_ATTACH
             if (InitCVIRTE (hinstDLL, 0, 0) == 0)
                             /* out of memory */
                 return 0;
         case DLL_PROCESS_DETACH:
             CloseCVIRTE ();
             break;
                                                                                                                                       |Waveform Graph|
     return 1:
  DLLEXPORT int acquisition(int carte, int voie, double* mesure)
  int status=0;
  //*mesure=8.3
  status = AI_VRead (carte, voie, 1, mesure);
                                                                                             agine on PCI-7030
  if(status!=0)
     return -1:
  return 0;
   39/41
              25 C
                      Ins I
```

可以使用C++寫好程式(或演算法),再使用LabVIEW呼叫

LabVIEW優點列舉(2/6):範例程式




- ■數千個範例程式
- ■LabVIEW內含的眾 多範例
- ■與網路整合的搜尋 工具

LabVIEW優點列舉(3/6):技術說明文件


- ■函式定義
- ■説明文件






LabVIEW優點列舉(4/6):遠端面板

網路



LabVIEW 環境




網頁瀏覽器

LabVIEW優點列舉(5/6): <u>分析資料</u>

內建超過500個以上訊號處理、訊號分析及數學運算的函式功能:

- ■訊號合成 (Signal synthesis)
- ■曲線近似 (Curve fitting) 及內插 (interpolation)
- ■利用 FFT 進行頻率分析
- ■數學運算
- ■機率與統計運算
- ■時域及頻域分析
- ■數位訊號處理
- ■波形排列及重新取樣
- ■更多其他的功能...
- ■針對量測的分析函式庫
- ■應用導向的工具箱




LabVIEW優點列舉(6/6):利用 LabVIEW 呈現

資料

資料視覺化

- 內建的使用者介面開發物件
- 資料趨勢圖及分佈圖
- 遠端應用監控





- 製作文件工具
- HTML 格式報表,可用網站瀏覽
- Microsoft Word 及 Excel 格式報表
- 與 NI DIAdem 軟體搭配產生互動式的 報表

資料管理及連接

- NI DIAdem 軟體
- 資料庫連結
- 存檔讀檔
- NI DataSocket 網路功能
- 開放的軟體環境 (ActiveX、DLLs)
- 標準的資料溝通方式
 - TCP/IP
 - □ 藍芽
 - OPC
 - □ SQL 資料庫連結
 - □ XML 資料格式

使用LabVIEW的優勢



- ■針對工程師及研究人員,設計的互動式圖形化程式語言
- ■提供數百個內建的函式,可進行資料輸入輸出、控制、分析及呈現等功能
- ■提供專門應用的高階開發工具及函式庫
- ■可將開發好的程式載入桌上型電腦、攜帶型電腦、工業電腦 及嵌入式處理器執行
- ■LabVIEW幾乎可以和所有的儀器進行溝通

課程目標

這個課程將可提供:

- ■了解人機介面、程式區、連接器以及常用的icon用法
- ■學習使用迴圈、以及常見的資料型別
- ■熟悉編輯與除錯技巧
- ■自行設計一個VI
- ■以圖表顯示資料以及儲存資料
- ■標準程式設計流程
- ■建構邏輯化的人機介面
- ■學習使用區域變數與廣域變數
- ■遠端資料轉裡

本課程不會提到以下幾點:

- ■DLL的使用或LabVIEW與C++的結合技巧
- ■設計儀器驅動程式
- ■GPIB或DAQ的詳細原理
- ■A/D, D/A原理

本課程需要使用的設備



執行LabVIEW8.2的個人電腦 建議使用Windows XP



